



Flask, Erlenmeyer, polymethylpentene capacity 125 mL

Številka izdelka:	2089843
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specifikacije

Enota:	1 each
Kapaciteta:	125 mL
Material:	Polymethylpentene
Vrsta:	Erlenmeyer