



Be Right™



Flask, Erlenmeyer, polymethylpentene capacity 250 mL, 4/pk

Številka izdelka: 2089876

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specifikacije

Enota:	4/pk
Kapaciteta:	250 mL
Material:	Polymethylpentene
Vrsta:	Erlenmeyer