



Pipeta, serološka, sterilna, polistiren, 5 mL, modra, 200 kosov

Številka izdelka:	2092667
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Sterilna plastična pipeta za enkratno uporabo. Posamično pakiranje. Barvno kodiranje. Umerjanje na količino TD. Skladne s standardom ASTM E934.

Graduations are calibrated for accurate dispensing to within $\pm 2\%$

Pipets meet BST/TSE requirements; latex-free

DNase/RNase-free, non-pyrogenic, non-cytotoxic, non-haemolytic

Specifikacije

Graduations:	0.10 mL
Kapaciteta:	5 mL
Negative Graduations:	2.0 mL
Prostornina / velikost paketa:	200 /pk