



Pipet rack

Številka izdelka:	2094700
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Polypropylene, fifty-place rack for safe, orderly storage of pipets. Helps create an efficient storage system or a well organized liquid handling station. Holes in top and middle decks are 16 mm diameter (0.63"); bottom holes are tapered to 6 mm (0.24") to prevent contact with supporting surface. Autoclavable to 122°C. Overall dimension: 21 x 11 x 22 cm (8.3 x 4.3 x 8.7").

Specifikacije

Dolžina:	210 mm
Globina:	22 cm
Prostornina / velikost paketa:	1 each
Širina:	4.3 Inches