



Crucible, Gooch, 40 mL Capacity

Številka izdelka: 2098641

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Specifikacije

Kapaciteta: 40 mL

Prostornina / velikost paketa: 1 /pk