



Reagent powder pillows, CuVer 1 copper, 10 mL, pk/1000

Številka izdelka:

2105828

Nearno

EUR Cena:

Kontakt

Na voljo



For copper determination by the bichinchoninate method. Range: up to 5.00 mg/L Cu. EPA compliant for wastewater reporting - Hach Method 8506. Sample Size: 10 mL. Pack of 1000 powder pillows.

Specifikacije

EPA compliant:	Da
Ime metode:	Bichinchoninate
Inštrument:	All except DR/3000, DR/2000, DR/850, DR/820
Merilno območje:	0.04 - 5.00 mg/L Cu
Parameter:	Copper
Platforma:	Praškasti reagenti
Prostornina / velikost paketa:	1000 /pk
Številka metode:	8506
Število testov:	1000