



Buffer, Phosphate Type, pH 7, for Amperometric Chlorine, 100 mL MDB

Številka izdelka: 2155332

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

For determination of chlorine by amperometric titration. APHA formulation 100 mL bottle

Specifikacije

Ime metode:	Amperometric Titration
Inštrument:	AutoCat [™] 9000
Merilno območje:	0,0051 - 5,000 mg/L Cl ₂
Platforma:	Raztopina
Številka metode:	8334
Število testov:	100