



Test kit, nitrite, nitrogen, model NI-15

Številka izdelka: 2182000

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Simple, low-cost color disc method is the most accurate visual test kit method.

Specifikacije

Ime metode:	Color disc/Diazotization
Merilno območje:	Single Range
Model:	NI-11
Parameter:	Nitrite
Platforma:	Color Disc
Smallest increment:	0.01
Standardno merilno območje:	0 - 0.40 mg/L NO ₂ ⁻ -N
Število testov:	100
Test Kit model:	NI-15
Vrsta primera:	D