



Test kit, formaldehyde, model FM-1

Številka izdelka:

2183100

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilises drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Drop count titration/ Thymolphthalein
Merilno območje:	0,05 - 1% & 0,5 - 10%
Parameter:	Formaldehyde - As CH ₂ O
Platforma:	Drop Count Titration
Smallest increment:	0.05% 0.5%
Število testov:	100
Test Kit model:	FM-1
Vrsta primera:	D