



Calcium Hardness Standard Solution, 10000 mg/L CaCO₃ (NIST), 16 x 10 mL Ampules

Številka izdelka:
EUR Cena:
Odprema v 2 tednih

218710
Kontakt

Nearno



For accuracy checks in calcium hardness determinations.

Simple accuracy checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the Tensette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimise waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specifikacije

Koncentracija:	10000 mg/L
Opis:	Calcium Standard Solution
Platforma:	Raztopina
Prostornina / velikost paketa:	10 mL /