



Filter, TSS glass fiber, pore size 1.5 µm, diameter 30 mm, 100/pk

Številka izdelka: 2198730

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Specifikacije

Diameter:	30 mm
Material:	Glass Microfiber
Prostornina / velikost paketa:	100 /pk
Velikost por:	1.5 µm