



Test kit, general boiler water, drop count titration, model AL-94 MG-L

Številka izdelka: 222403
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Used for controlling boiler and cooling water treatment and checking corrosion, this test kit is an economical solution for the small boiler operator who needs to monitor alkalinity, chloride, hardness, iron, and pH.

Uses simple, economical drop count titration and color disc methods

Kit contains all required reagents and apparatus to perform analysis in the field

Uses drop count titration and color disc methods

Specifikacije

Aplikacija:	Boiler Water
Model:	AL-94 MG-L
Parameter:	Drop Count Titration
Platforma:	Color Disc