



Bag, sterile, WhirlPak, 207 mL, 25/pk

Številka izdelka:	2233125
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

Specifikacije

Kapaciteta:	207 mL
Prostornina / velikost paketa:	25