



Test kit, nitrogen , model NI-8

Številka izdelka: 224100
EUR Cena: Kontakt
Datum pošiljanja ni prikazan

Nearno



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikacije

Ime metode:	Color disc/Nessler reagent
Merilno območje:	0 - 3.0 mg/L NH ₃ -N
Parameter:	Nitrogen, Ammonia, mid range - For freshwater
Platforma:	Color Disc
Smallest increment:	0.1
Število testov:	100
Test Kit model:	NI-8
Vrsta primera:	B