



## Test kit, chlorine (total), model CN-65

Številka izdelka:

225401

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifikacije

Ime metode:	Drop count titration/ Thiosulfate
Merilno območje:	0,2 - 4,0 mg/L & 1 - 20 mg/L
Najmanjši korak rasti:	1
Parameter:	Chlorine, mid and high range - As Total Cl <sub>2</sub>
Platforma:	Drop Count Titration
Smallest increment:	0.2
	1
Število testov:	100
Test Kit model:	CN-65
Vrsta primera:	D