



## Sulfite Equivalent Standard Solution, 5000 mg/L as SO<sub>3</sub>, pk/16 - 10 mL Voluette Ampules

Številka izdelka: 2267410

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Sodium Thiosulfate Standard Solution (0.125 N), equivalent to 5000 mg/L as SO<sub>3</sub>

### Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

---

## Specifikacije

|                                |   |
|--------------------------------|---|
| Koncentracija:                 | 5000 mg/L   |
| Opis:                          | Sulfite Equivalent Solution as SO <sub>3</sub> <sup>2-</sup> (NIST) |
| Platforma:                     | Raztopina   |
| Prostornina / velikost paketa: | mL/16   |