



Buffer Solution, pH 7.00, Colour-coded Yellow, 20 L

Številka izdelka: 2283561
 EUR Cena: Kontakt
 Odprema v 2 tednih

pH Buffer standard solutions for laboratory and field analyses

pH Buffer Standard Solution for pH electrode calibration. Colour-coded pH buffer solution for easy identification. pH is 7.00 ± 0.02 at 25°C . 20 L Bag-in-box package with a collapsible inner bladder in a white corrugated cardboard box. Ideal packaging size for laboratories doing high volume pH testing.

Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Barva kodirana?:	Da
Koncentracija:	$7,00 \pm 0,02$ pH
Opis:	pH Buffer Calibration reagent
Package settings:	Bag-in-box (BiB)
pH vrednost(i):	7.00 ± 0.02 pH at 25°C
Pogoji skladiščenja:	$10 - 25^{\circ}\text{C}$ (protect from light)
Posebni zapiski:	Observed shelf-life after opening, following recommended storage and use: 6 months
Prostornina / velikost paketa:	20 L
Sledljivo po NIST-u:	Da