



## Buffer Powder Pillows, pH 4.01 & pH 7.00, pk/2x15

Številka izdelka: 2299264  
EUR Cena: Kontakt  
Datum pošiljanja ni prikazan

Nevarno



### pH Buffer standard solutions for laboratory and field analyses

pH Buffer Powder Pillows for pH electrode calibration. Each powder pillow prepares 50 mL of fresh buffer solution. Colour-coded for easy identification. pH is  $\pm 0.02$  at 25°C. Pack of 15 pillows of pH 4.01 (red) and 15 pillows of pH 7.00 (yellow).

#### Traceable to NIST standard reference materials

Whether in the laboratory or in the field, always perform measurement results you can trust is essential. The quality of the standard pH Buffer Solutions that are used are key.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

#### Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### Specifikacije

|                                |                                     |
|--------------------------------|-------------------------------------|
| Barva kodirana?:               | Da                                  |
| Koncentracija:                 | 4,01 pH and 7,00 pH                 |
| Opis:                          | pH Buffer Calibration reagent       |
| Package settings:              | Powder Pillows                      |
| pH vrednost(i):                | 4.01 & 7.00 $\pm 0.02$ pH at 25°C   |
| Pogoji skladiščenja:           | 10 - 25 °C (protect from light)     |
| Posebni zapiski:               | Discard prepared solution after use |
| Prostornina / velikost paketa: | 30 pcs                              |
| Sledljivo po NIST-u:           | Da                                  |