



Test kit, iron color disc, model IR-21

Številka izdelka:

2299300

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



Iron color disc/TPTZ low range test kit contains 100 tests.

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/TPTZ
Merilno območje:	0 - 2.0 mg/L Fe
Model:	IR-21
Parameter:	Iron TPTZ
Platforma:	Color Disc
Smallest increment:	0.01
	0.05
Standardno merilno območje:	0 - 2.0 mg/L Fe; 0 - 1.2 mg/L Fe
Število testov:	100
Test Kit model:	IR-21
Vrsta:	Double Metric
Vrsta primera:	D