



Standard solution, sodium chloride, 85.47 mg/L (180 µS/cm), 1L

Številka izdelka: 2307553

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Conductivity standard solutions for laboratory and field analyses

For standardization and accuracy checks in conductivity measurements

Specifikacije

Formula:	NaCl
Koncentracija:	180 µS/cm at 25°C ±10 µS/cm at 25°C
Opis:	Conductivity calibration solution
Package settings:	Bottle
Parameter:	Conductivity
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	1000 mL