



Filter, acrocap, pore size 0.2 µm

Številka izdelka: 2323810

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

Specifikacije

Diameter:	NA
Filter Type:	Filter, acrocap
Packing:	NA
Velikost por:	0.2 µm