



Bottle, sampling, polycarbonate, 500 mL, 4/pk

Številka izdelka:	2324539
EUR Cena:	Kontakt
Datum pošiljanja ni prikazan	

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifikacije

Body Material:	Polycarbonate
Kapaciteta:	500 mL
Prostornina / velikost paketa:	4 /pk