



Be Right™



Bottle, Sampling, Polycarbonate, 500 mL, 24/pk

Številka izdelka: 2324579

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifikacije

Body Material: Polycarbonate

Kapaciteta: 500 mL

Prostornina / velikost paketa: 24 /pk