



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

Številka izdelka: 2324663

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifikacije

Body Material:	Polycarbonate
Kapaciteta:	1000 mL
Prostornina / velikost paketa:	4 /pk