



Fluoride Standard Solution, 100 mg/L, 500 mL

Številka izdelka:	23249
EUR Cena:	Kontakt
Na voljo	

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Fluoride, 100 mg/L ± 2 mg/L F⁻ (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Koncentracija:	100 mg/L ± 2 mg/L F ⁻
Opis:	Fluoride Standard Solution
Package settings:	Bottle
Parameter:	Fluoride
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL