



Filter holder, for vacuum filtration, glass

Številka izdelka: 234000

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

Specifikacije

Diameter:	47 mm
Prostornina / velikost paketa:	1
Vsebina paketa:	Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).

Vsebina paketa

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).