



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Številka izdelka: 2349500

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Nearno



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikacije

Ime metode:	Color disc/Iron reduction
Merilno območje:	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Model:	OS-1
Parameter:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platforma:	Color Disc
Število testov:	50
Test Kit model:	OS-1
Vrsta primera:	R