



## Standard solution, potassium, 100 mg/L as K (NIST), 500 mL

Številka izdelka: 2351749

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

### **Ion Selective Electrode (ISE) standard solutions for laboratory analyses**

For accuracy checks in potassium determinations

#### **Traceable to NIST standard reference materials**

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

#### **Design for use with Hach meters portfolio and other meters.**

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

#### **Reliable packaging**

Supplied in an airtight packaging ensuring guaranteed shelf-life.

---

### **Specifikacije**

Koncentracija:	100 mg/L $\pm$ 5 mg/L K <sup>+</sup>
Opis:	Potassium Standard Solution
Package settings:	Bottle
Parameter:	Potassium
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	500 mL