



Ammonia Standard Solution, 1000 mg/L, 1 L

Številka izdelka: 2354153
EUR Cena: Kontakt
Na voljo

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Ammonia, 1000 mg/L ± 20 mg/L $\text{NH}_3\text{-N}$. 1 L bottle

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifikacije

Koncentracija:	1000 mg/L ± 20 mg/L $\text{NH}_3\text{-N}$
Opis:	Ammonia Standard Solution
Package settings:	Bottle
Parameter:	Ammonium
Platforma:	Raztopina
Pogoji skladiščenja:	10 - 25 °C (protect from light)
Posebni zapiski:	Observed shelf life after opening, following recommended storage and use: 3 months
Prostornina / velikost paketa:	1 L