



## Test kit, nitrogen, ammonia, model NI-SA

Številka izdelka:

2428700

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

**Nearno**



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

**Designed for testing ammonia in freshwater or seawater**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifikacije

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Ime metode:         | Color disc/Salicylate                                       |
| Merilno območje:    | 0 - 2.4 mg/L NH <sub>3</sub> -N                             |
| Parameter:          | Nitrogen, Ammonia, mid range - For freshwater or seawater** |
| Platforma:          | Color Disc                                                  |
| Smallest increment: | 0.1                                                         |
| Število testov:     | 100                                                         |
| Test Kit model:     | NI - SA                                                     |
| Vrsta primera:      | D                                                           |