



Reagent set, hardness (total), HA-71A test kit

Številka izdelka:

2437800

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For determination of Total Hardness by the Drop Count EDTA Titration Method. Replacement reagent set for HA-71A Test Kit. Set contains Buffer Solution (Hardness 1), Hardness 2 and Hardness 3 Test Solutions for approximately 100 tests.

Specifikacije

Ime metode:	EDTA Titration
Merilno območje:	1 - 20 mg/L CaCO ₃
Parameter:	Hardness (Total)
Posebna značilnost:	Reagent set, HA-71A test kit