



## Reagent set, copper, free & total

Številka izdelka: 2439200  
EUR Cena: Kontakt  
Datum pošiljanja ni prikazan

Nearno



For determination of free and complexed (chelated) Copper by the Binchoninate method. Range: up to 5.00 mg/L. Set includes: Free Copper Reagent Powder Pillows and Hydrosulfite Power Pillows for 100 tests.

### Specifikacije

|                  |                                         |
|------------------|-----------------------------------------|
| Ime metode:      | Bicinchoninate                          |
| Inštrument:      | DR 6000, DR 3900, DR 1900, DR 900, PCII |
| Merilno območje: | 0.04 - 5.00 mg/L Cu                     |
| Platforma:       | Praškasti reagenti                      |
| Število testov:  | 100                                     |