



Reagent set, chloride, LR, digital titrator

Številka izdelka:

2444100

EUR Cena:

Kontakt

Datum pošiljanja ni prikazan

Nearno



For determination of Chloride by argentometric Digital Titrator titration. Range: 10 - 100 mg/L Cl⁻. Reagent Set contains Chloride 2 Indicator Powder Pillows and Silver Nitrate Digital Titrator Cartridges (2) 0.2256 N. Approximately 100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikacije

Ime metode:	Silver Nitrate Titration
Merilno območje:	10 - 100 mg/L Cl ⁻
Številka metode:	8207
Število testov:	100