



Reagent set, hardness (total), buret, LR

Številka izdelka:

2447600

Nearno

EUR Cena:

Kontakt

Odprema v 2 tednih



Reagent set for the determination of Total Hardness by the ManVer[®] 2 Buret Titration Method (Hach method 8226). EPA Accepted method. Range: 0-2,500 mg/L as CaCO₃. Contains: Buffer Solution (Hardness 1), ManVer[®] 2 Hardness Indicator Powder Pillows and 0.0200 N EDTA (TitraVer[®]) Titrant for up to approximately 100 tests.

Specifikacije

EPA compliant:	Da
Ime metode:	EDTA Titration
Merilno območje:	0 - 2500 mg/L CaCO ₃
Parameter:	Hardness (Total)
Posebna značilnost:	Reagent set, Buret, LR
Številka metode:	8226
Število testov:	100