



Reagent set, hardness (total), buret, HR

Številka izdelka:

2447700

Nearno

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.



Reagent set for the determination of Total Hardness by the ManVer® 2 Buret Titration Method. Range: 100-25,000 mg/L as CaCO₃. Contains: Buffer Solution (Hardness 1), ManVer® 2 Hardness Indicator Powder Pillows and 0.200 N EDTA (TitraVer®) Titrant for up to approximately 100 tests.

Specifikacije

EPA compliant:	Da
Ime metode:	EDTA Titration
Merilno območje:	100 - 25000 mg/L CaCO ₃
Parameter:	Hardness (Total)
Posebna značilnost:	Reagent set, Buret, HR
Številka metode:	8222
Število testov:	100