



Flask, Erlenmeyer, Wide Mouth, 2000 mL, Subdivisions 600 mL

Številka izdelka: 2489454

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included. Requires Size 13 rubber/silicone stopper.

Specifikacije

Enota:	1 each
Kapaciteta:	2000 mL
Material:	Glass Wide Mouth