



Bottle, sampling, presterilized, round, 12/pk

Številka izdelka: 2495012

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Works well with P-A Broth Ampules and PathoScreen Medium.

Specifikacije

Dechlorinating Agent?:	No
Prostornina / velikost paketa:	12/pk