



Sample Cell: 1" Square Glass 10 mL matched set of 8

Številka izdelka: 2495408

EUR Cena: Kontakt

Na voljo

Comes as a matched set of 8 cells.

Use in Combination with Hach Spectrophotometers

Hach spectrophotometers, chemistries, and methods are designed, developed, and manufactured to optimize your results.

Specifikacije

Komentar:	1 pair of cells is included with purchase of DR 5000 UV-Vis, DR 3800, DR 2800, or DR 2700 Spectrophotometer
Prostornina / velikost paketa:	8 /pk