



Rack, Test Tube, Polyethylene, 13 mm (O.D.), 90 Holes

Številka izdelka: 2497900

EUR Cena: Kontakt

Odprema v 2 tednih

Ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specifikacije

Asortiman testov:	6 x 15
Diameter:	13 mm
Luknje za teste:	90
Prostornina / velikost paketa:	1 each