



Test tube, rack, polyethylene, 16 mm (O.D.), 60 holes

Številka izdelka: 2497901

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Sixty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may vary.

Specifikacije

Asortiman testov:	5 x 12
Diameter:	16 mm
Luknje za teste:	60
Prostornina / velikost paketa:	1 each