



Test Tube, rack, polyethylene, 20 mm (O.D.), 20 holes

Številka izdelka:	2497912
EUR Cena:	Kontakt
Na voljo	

Twenty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may Vary.

Specifikacije

Asortiman testov:	4 x 5
Diameter:	20 mm
Luknje za teste:	20
Prostornina / velikost paketa:	1 each