



**Be Right™**



## ORP buffer solution, 600 mV, 3.5 L

Številka izdelka: 25M2A1002-123

EUR Cena: Kontakt

Odprema v 2 tednih

ORP Reference Solution used for ORP sensor calibration (in resealable plastic bottle).

---

### Specifikacije

|                                |                          |
|--------------------------------|--------------------------|
| Inštrument:                    | Digital Differential ORP |
| Merilno območje:               | 600 mV                   |
| Parameter:                     | ORP                      |
| Prostornina / velikost paketa: | 3.5 L                    |