



Filter, membrane, pore size 0.45 µm, diameter 25 mm, 100/pk

Številka izdelka:

2514101

Nearno

ZASTAREL IZDELEK
Ta izdelek ni več na voljo.



For use in a wide-range of laboratory applications, these discs are composed of Mixed-Cellulose Esters and are the industry standard for particulate removal, particulate analysis and clarification of fluids. Diameter: 25mm. Pore Size: 0.45µm. 100/PK.

Specifikacije

Diameter:	25 mm
Velikost por:	0.45 µm