



Be Right™



Test Kit, AccuVac® Color Disc, Ozone, 0.05-1.50 mg/L

Številka izdelka:
EUR Cena:
Datum pošiljanja ni prikazan

2518050

Kontakt

Nevarno



This ozone test kit uses the fast, accurate visual method for determining ozone concentrations. Use to control disinfectant levels in water, wastewater, and industrial water. Kit contains 25 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

AccuVac® vials mean no mess, no measuring

Indigo trisulfonate method masks interferences from chlorine and other oxidizing agents and makes ozone analysis easier

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Specifikacije

Ime metode:	Color disc/ Indigo trisulfonate
Merilno območje:	0 - 1.50 mg/L O ₃
Model:	NA
Parameter:	Ozone
Platforma:	Accuvac
Smallest increment:	0.05
Standardno merilno območje:	0 - 1.50 mg/L O ₃
Število testov:	25
Test Kit model:	AccuVac® Kit
Vrsta primera:	G