



Bottle rack, plastic-coated, 20 compartments

Številka izdelka: 2551801

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak® sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Specifikacije

Posebna značilnost:	Rack, 20 Bottle Compartments
Prostornina / velikost paketa:	1