



Sample container, 240 mL

Številka izdelka: 2552411

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Plastic, leak-proof. Corning Snap-Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with unique sealing mechanism so sample won't leak.

Specifikacije

Kapaciteta:	240 mL
Premer, zunanji:	80 mm (approx)
Prostornina / velikost paketa:	140 /pk