



Rack, test tube, no-wire

Številka izdelka: 2587300

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

Ninety-place No-Wire Rack features locking mechanism that prevents “floaters” in water baths and speeds decanting by allowing the whole rack to be inverted. Heavier-than-water polypropylene assures non-floating racks when immersed, even in agitated water. Autoclavable. Numbered and lettered rows (6 x 15 array) aid in identification of tubes (10-13 mm diameter). Collapsible for storage, yet assemble in seconds. Overall dimensions: 24.8 x 10.5 x 7 cm (9.75 x 4.12 x 2.75").

Specifikacije

Dolžina:	248 mm
Globina:	7 cm
Prostornina / velikost paketa:	1 each
Širina:	4.12 Inches