



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Številka izdelka: 2618800

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Specifikacije

Diameter:	47 mm
Filter Type:	Filter
Material:	Polyethersulfone
Prostornina / velikost paketa:	200 /pk
Velikost por:	0.2 µm