



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Številka izdelka:	2634001
EUR Cena:	Kontakt
Odprema v 2 tednih	

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. (**Please Note:** These sample vials are not autoclaveable)

Specifikacije

Diameter:	43mm
Dolžina:	93mm
Kapaciteta:	fill line at 100 mL (+/-2.5 mL)
Posebna značilnost:	Sample Vial, with Thiosulfate
Prostornina / velikost paketa:	200 /pk
Unit:	200 /pk