



Be Right™



Sample vial, 120 mL w/out thiosulfate, 200/pk

Številka izdelka: 2634002

EUR Cena: Kontakt

Datum pošiljanja ni prikazan

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene.

Specifikacije

Diameter:	43 mm
Dolžina:	93 mm
Kapaciteta:	Fill line at 100 mL ($\pm$ 2.5 mL)
Posebna značilnost:	Sample Vial, without Thiosulfate
Prostornina / velikost paketa:	200/pk